FYZIOPONK

DESATERO

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FYZIOPONK

TEN COMMANDMENTS

- 1.** The purpose of using the rehabilitation pad Fyzioponk is to treat clinical conditions indicated by the physiotherapist or the attending physician for rehabilitation using the Vojta method and to treat and alleviate various deviations in psychomotor development. The indications for this rehabilitation are diagnoses such as cerebral palsy, torticollis, brachial plexus palsy, spina bifida, pes equinovarus, scoliosis, orthopaedic defects of the chest and other neurological diseases requiring intensive rehabilitation care of users of childhood (newborns, infants, toddlers and younger children).
- 2.** Always place the mat on a flat and clean surface (table, chest of drawers, kitchen counter, etc.).
- 3.** Never leave your child unattended on the mat! A fall and subsequent injury may occur.
- 4.** Always make sure the pad is on a stable surface and cannot tip over.
- 5.** The surface of the pad is washable. Before and after use, you can wipe it with soapy water or common detergents using a sponge or cloth. Cleaning and disinfection MAY be done with 98% alcohol. Do not use other than recommended cleaning agents - this could damage the material.
- 6.** Check the pad regularly to check for loose clips or damage to the leatherette or other parts.
- 7.** Transport the mat sufficiently packed. This will ensure a long service life.
- 8.** Keep the mat in a dry place and do not leave heavy objects on it.
- 9.** When exercising, always follow the advice and instructions of your physiotherapist or your doctor.
- 10.** Remember that the practice of the Vojta method is a process. Even the best mat will not do it for you and the child. Be patient.

The shelf life of the pad is at least 3 years. If the pad is damaged, do not use it. After the expiry date or in case of damage, dispose of the pad in municipal waste.

October 28, 2021